

Syndromic Surveillance in veterinary medicine using laboratory submission data – lessons learned from two systems

Fernanda Dórea

Ann Lindberg - SVA - Sweden

Beverly McEwen - AHL, UoG - Canada

Crawford Revie - UPEI - Canada

Javier Sanchez - UPEI - Canada

Information about the systems

- **Canada**
 - Animal Health Laboratory (AHL) in the province of Ontario
 - Laboratory of choice for livestock diagnostics
 - **Laboratory submissions (test orders) in cattle**
- **Sweden**
 - National Veterinary Institute
 - All testing of reportable diseases in animal in Sweden
 - **All laboratory submissions**

Records classification




caseID	Source	Date	Sample	Test
2002814	clinical	2000-jan-01	Serum	BLV
2002814	clinical	2000-jan-01	Serum	BLV
2002582	clinical	2000-jan-01	Serum	BVD type 2
2002607	clinical	2000-jan-01	Serum	BVD type 2
2002607	clinical	2000-jan-01	Serum	BVD type 2
2000032	clinical	2000-jan-02	Liver	Bacterial culture
2000032	clinical	2000-jan-02	Lung	Bacterial culture
2000032	clinical	2000-jan-02	Spleen	Bacterial culture
2000032	clinical	2000-jan-02	Tissue,_Fixed	Histopathology
2000128	clinical	2000-jan-02	Tissue,_Fixed	Histopathology
2000055	clinical	2000-jan-02	Tissue,_Fixed	BVD type 1
2000845	clinical	2000-jan-02	Brain	Bacterial culture
2000001	clinical	2000-jan-02	Brain	Bacterial culture
2015705	clinical	2000-jan-02	Fetal_Tissue	Bacterial Culture
2015789	clinical	2000-jan-02	bone	Clostridia
2015789	clinical	2000-jan-02	joint	Clostridia
2015719	clinical	2000-jan-02	Lung	Bacterial culture
2015705	clinical	2000-jan-02	rumen	Bacterial culture
2015705	clinical	2000-jan-02	rumen	Bacterial culture
2015705	clinical	2000-jan-02	Bronchoalveolar_Lavage	Bacterial culture
2024811	clinical	2000-jan-02	large int	Clostridium perfringens


Classification of records into syndromes

Syndrome	caseID	Source	Date	Sample	Test
BLV	2002814	clinical	2000-jan-01	Serum	BLV
BLV	2002814	clinical	2000-jan-01	Serum	BLV
BVDV	2002582	clinical	2000-jan-01	Serum	BVD type 2
BVDV	2002607	clinical	2000-jan-01	Serum	BVD type 2
BVDV	2002607	clinical	2000-jan-01	Serum	BVD type 2
Hepatic	2000032	clinical	2000-jan-02	Liver	Bacterial culture
Respiratory	2000032	clinical	2000-jan-02	Lung	Bacterial culture
Haematopoietic	2000032	clinical	2000-jan-02	Spleen	Bacterial culture
Non-specific	2000128	clinical	2000-jan-02	Tissue,_Fixed	Histopathology
Non-specific	2000128	clinical	2000-jan-02	Tissue,_Fixed	Histopathology
BVD	2000055	clinical	2000-jan-02	Tissue,_Fixed	BVD type 1
Nervous	2000845	clinical	2000-jan-02	Brain	Bacterial culture
Nervous	2000001	clinical	2000-jan-02	Brain	Bacterial culture
Abortion	2015705	clinical	2000-jan-02	Fetal_Tissue	Bacterial Culture
Musculoskeletal	2015789	clinical	2000-jan-02	bone	Clostridia
Musculoskeletal	2015789	clinical	2000-jan-02	joint	Clostridia
Respiratory	2015719	clinical	2000-jan-02	Lung	Bacterial culture
GIT	2015705	clinical	2000-jan-02	rumen	Bacterial culture
GIT	2015705	clinical	2000-jan-02	rumen	Bacterial culture
Respiratory	2015705	clinical	2000-jan-02	Bronchoalveolar_Lavage	Bacterial culture
GIT	2024811	clinical	2000-jan-02	large int	Clostridium perfringens
Respiratory	2024899	clinical	2000-jan-02	Lung	Virus isolation


Case counts for each syndrome




Retrospective analysis



Prospective analysis



Alarm detection + Real-time monitoring



1. Classification: records → syndrome counts

- *Lack of standardized coding systems*
- Challenges:
 - Limited clinical information
 - *Active versus passive activities*
 - Structured limited information + free-text
 - Missing or erroneous information

1. Classification: records → syndrome counts

- *Exhaustive syndromic list*
- Methods
 - Direct mapping
 - Rule-based - 99.4% recall and 100% precision
 - Dictionary of relevant words
 - Direct relationships
 - Co-occurrences
 - Precedence

Animal groups

Cattle

Equine

Swine

Poultry

Cats

Dogs

Aquatic

Wild

Environmental

Feed

Syndromic groups

Abortion

Circulatory, Hepatic, Haematopoietic

Eyes, Ears and Others

Integumentary

GIT

MBS

Mastitis

Musculoskeletal

Nervous

Reproductive

Respiratory


Systemic

Urinary

AntimResistance

Fallen stock

2. Retrospective analysis




- Noise versus past outbreaks
- Temporal effects
- *Small counts per day*

2. Retrospective analysis


AHL

Systemic




SVA

Abortion




SVA

Abortion



3. Prospective analysis



3. Prospective analysis

- What algorithms
 - Control charts (Shewhart + EWMA)
 - *R:{qcc}*
 - Pre-processing: Differencing or Regression Models
 - Holt-Winters exponential smoothing
 - *R:HoltWinters()*


Characteristics
of the
outbreak

Characteristics
of the data
stream


Respiratory - 2011-12-09

Algorithm	D-4	D-3	D-2	D-1	TODAY
HW	0	0	0	0	5
EWMA	0	0	0	0	2
Shewhart	0	3	0	0	5
SUM	0	3	0	0	12












Respiratory DAILY Alarm Score



Respiratory Observed counts and Corrected Baseline



[Main page](#)

-  [Cats](#)
-  [Cattle](#)
-  [Dogs](#)
-  [Environm.](#)
-  [FOD](#)
-  [FSK](#)
-  [Horses](#)
-  [Poultry](#)
-  [Small Ruminants](#)
-  [Swine](#)
-  [VLT](#)

Syndromic surveillance at SVA

2013-07-27

Select a group on the navigation menu to the left to see results

System outputs are based on data up to the end of the PREVIOUS day.

Number of events per day correspond to the number of laboratory submissions, classified into syndromic groups by a computer system.

They reflect the number of ANIMALS tested per day for CATS and DOGS, and the number of HERDS otherwise


For questions contact Fernanda Dorea (fernanda.dorea@sva.se)

Free tools




Basic framework

SYNDROMIC SURVEILLANCE




Enriching the system INFORMATION SYSTEM



Acknowledgments



STATENS
VETERINÄRMEDICINSKA
ANSTALT



fernanda.dorea@sva.se

"IN GOD WE TRUST;
ALL OTHERS MUST
BRING DATA."
- W. EDWARDS DEMING